



MelMedtronics Holdings, LLC

David W. Holmes, Ph.D.

CEO & Senior Scientist

Board Certified in Audiology


2083 Rana Park
Flint, TX 75762
817-201-6791 Cell
903-705-7450 Fax
drdavidholmes Skype

holmes@melmedtronics.com

www.melmedtronics.com

TINNITUS PREVALENCE

600M+ People Worldwide



- ▶ 15% of world-wide population has has tinnitus
- ▶ 26.7% people ages 65-84 years suffer tinnitus
- ▶ 55M just in USA

TINNITUS: ISSUES

Clinicians who see patients with tinnitus share a common axiom:

“If you’ve seen one tinnitus patient, you’ve seen one tinnitus patient.”

Patients with tinnitus are very complex for more reasons than can possibly be listed, but here are a few examples:

- 600+ prescription drugs have a tinnitus side effect
- Tinnitus effects of interactions among multiple drugs are unknown
- 200+ known medical conditions
- Variety of frequency and intensity parameters of noise induced tinnitus

TINNITUS: VARIABLES

- **Personality Type Reactions (Jung)**
 - Sensing-Judging
 - Sensing-Perceiving
 - Intuitive-Feelers
 - Intuitive-Thinkers
- **Anxiety**
- **Depression**
- **Hyperacusis:** reduction of normal tolerance for everyday sounds
- **Misophonia:** hatred of sound
- **Phonophobia:** pathological fear of sound, or of speaking aloud

TINNITUS: CURRENT SOLUTIONS

CURRENT TREATMENTS

- Sound therapy: Maskers
- Hearing Aids: amplify ambient noise
- Homeopathic

SUCCESS?

- No cures
- Non-compliant patients
- Marginal success for some patients

ISSUES

- Patients' tinnitus may change and fluctuate over the time
- If the treatment is partially effective, it may change the perception of the frequency and/or the intensity of the patient's tinnitus for a variety of reason

TINNITUS NEUROMODULATION TREATMENT

TNT™

- ▶ Develop a **“modifiable” treatment program**, that would respond to changes in the patient’s tinnitus
- ▶ The use of **“intelligent” algorithms** could then be used to predict which program had the highest probability of success, based on variables that had been identified through ongoing diagnostic information

PATIENT DATA

- **Demographics**
- **History**
 - Date of onset
 - Etiology
 - Medications
 - Medical conditions

ASSESSMENTS

- **Description of Tinnitus**
 - **Types**
 - **Tonal**
 - Frequency
 - Minimum masking level
 - Residual inhibition
 - **Noise**
 - Intensity
 - Minimum masking level
 - Residual inhibition
 - **Location**
 - Left, Right, Center

TNT™: ASSESSMENTS

- **Psyco/Social Measures**
 - Tinnitus Handicap inventory
 - Tinnitus Reaction Quotient
 - Tinnitus Severity Scale
 - Depression scale
 - Anxiety scale
- **Personality Index (Jung)**
 - Sensing-Judging
 - Sensing-Perceiving
 - Intuitive-Feelers
 - Intuitive-Thinkers

These will be obtained using surveys

TNT™: ACOUSTIC TREATMENT OPTIONS

- **Tonal**

- Non-modulated frequency
- Frequency Modulated (40 Hz)
- Frequency Modulated (Sweep 40-100 Hz)
- Amplitude Modulated (+/- 5 dB)
- Combined Modulated (AM & 40 Hz)

- **Noise**

- Non-modulated Notched band
- Modulated notched band (40 Hz)
- Modulated notched band (sweep 40-100 Hz)
- Notched band (AM +/- 5 dB)
- Combined notched band (AM and FM 40 Hz)
- Combined notched band (AM and sweep FM 40-100 Hz)

- **Zwicker**

- Non-modulated notched Zwicker band
- Modulated notched Zwicker band (40 Hz)
- Modulated notched Zwicker band (sweep 40-100 Hz)

TNT™: ACOUSTIC TREATMENT OPTIONS

- **Music**

- Non-modulated Notched band
- Modulated notched band (40 Hz)
- Modulated notched band (sweep 40-100 Hz)
- Notched band (AM +/- 5 dB)
- Combined notched band (AM and 40 Hz)
- Combined notched band (AM and sweep 40-100 Hz)

- **Fractal**

- Non-modulated Notched band
- Modulated notched band (40 Hz)
- Modulated notched band (sweep 40-100 Hz)
- Notched band (AM +/- 5 dB)
- Combined notched band (AM and FM 40 Hz)
- Combined notched band (AM and sweep FM 40-100 Hz)

TNT™: PROCEDURES

PATIENT REGISTRATION

- Name
- Date of birth
- Sex
- Unique I.D.
- Login (email address)
- Password

DATA

- **Data logging**
 - Login/logout date, time and duration treatment sessions
 - Treatment assigned
- **Data collection**
- **Analyses**

TNT™: PROCEDURES PHASE I

- **Match tone frequency**
 - generate a narrowband of +/- 10% of matched tone
- **Match narrowband noise**
 - Modulate frequency of tonal band 40 Hz
 - Sweep from 40 Hz to 100 Hz back and forth over a 5 sec. window
- **Modulate amplitude**
 - +/- 5 dB
- **Zwicker - Filter out notched frequency**
 - (+/- 10% of matched tone)
 - Noise narrowband

SELECTED TREATMENTS

- Amplitude modulation +/- 5 dB of threshold
- Frequency modulation
 - 40 Hz – 100 Hz
 - 40 Hz
- Zwicker

TNT™: PROCEDURES PHASE II

CLINICAL TRIALS

- Determine N subjects based on number of variables
- Contact potential sites
- Recruit subjects through:
 - Social Media
 - Tinnitus groups
 - Website
- **Develop**
 - data collection procedures
- **Integrate**
 - Data
 - Statistical recommendations
- **Select**
 - Program selection



MelMedtronics Holdings, LLC

Where Ideas Become Reality

www.melmedtronics.com